HARTMANN

HARTMANN

Tensoval® mobil

Gebrauchsanleitung
Mode d'emploi
Gebruiksaanwijzing
Instructions for use
Instrucciones de uso
Manual de utilização
تعليمات الاستخدام

PAUL HARTMANN AG · 89522 Heidenheim, Germany

AE – PAUL HARTMANN Middle East FZE · Dubai

AT – PAUL HARTMANN Ges.m.b.H. · 2355 Wiener Neudorf

BE – N.V. PAUL HARTMANN S.A. · 1480 Saintes/Sint-Renelde

CH – IVF HARTMANN AG · 8212 Neuhausen

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Preliminary remarks

⚠ Please read these instructions carefully before first use as correct blood pressure measurement depends on the appropriate use of the device. These instructions for use are designed to instruct you, from the very start, in the individual steps of self-measurement of blood pressure using Tensoval mobil. You will thus receive important and helpful hints for producing reliable results for your personal blood pressure profile. Be sure to keep these instructions for use for future reference.

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NL – PAUL HARTMANN B.V.
Postbus 26
6500 AA Nijmegen

BE – N.V. PAUL HARTMANN S.A.
1480 Saintes/Sint-Renelde

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1. Introduction	Page 82
2. General information on blood pressure	82
2.1 WHO and ISH classification of high blood pressure limits	82
2.2 The importance of self-measurement of blood pressure	83
2.3 Objectives of self-measurement of blood pressure	84
2.4 Regular blood pressure measurement	85
3. Getting ready for self-measurement	85
3.1 Inserting the batteries	85
3.2 Setting date and time	86
3.3 Correct posture and position	87
3.4 Applying the wrist monitor	87
4. Measuring blood pressure	88
5. Setting the memory function	90
5.1 Saving the measured values	90
5.2 Recalling the measured values	90
5.3 Deleting the measured values	91
5.4 Using the guest mode	92
6. Maintenance of the device	92
7. Tensoval range of products	93
8. Warranty conditions	93
9. Explanation of error displays	94
10. Important notes	97
10.1 Drugs	97
10.2 Pregnancy	98
10.3 Diabetes, history of other medical conditions	98
10.4 Arrhythmias, heart rhythm disorders (cardiac arrhythmias), cardiac pacemakers	98

11. Control displays and symbols	Page 99
12. Technical data	100
13. Power supply, disposal note and safety information	101
13.1 Batteries, and disposal	101
13.2 Safety information	101
14. Legal requirements and guidelines	102
15. Measurement function check and service addresses	102
15.1 Explanation of measurement function check	102
15.2 Contact information for customer queries	103

1. Introduction

Dear Customer,
We are pleased that you have decided to purchase a wrist blood pressure monitor from HARTMANN. Tensoval mobil is a quality product for fully automatic self-measurement of blood pressure on the wrist. Requiring no preliminary settings, this device conveniently inflates automatically for easy, quick and reliable measurement of the systolic and diastolic blood pressures as well as the pulse rate. Using the HARTMANN Comfort Air Technology, this device automatically determines the maximum inflation pressure to be used for comfortable individualised blood pressure meas-

urements. This device is designed to optimally support you in taking blood pressure measurements. We wish you all the best for your health.

2. General information on blood pressure

2.1 WHO and ISH classification of high blood pressure limits

The World Health Organisation (WHO) and the International Society of Hypertension (ISH) have developed the following classification for blood pressure values:

World Health Organisation (WHO) blood pressure targets, 1999

Assessment	Systolic pressure	Diastolic pressure
Optimal	up to 120 mmHg	up to 80 mmHg
Normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg

To determine your blood pressure you need to measure two values:

- Systolic (upper) blood pressure: is produced when the heart contracts and expels blood into the blood vessels.
- Diastolic (lower) blood pressure: this is the value measured when the heart muscle is dilated and again fills with blood.
- Blood pressure readings are expressed in mmHg.

Established hypertension (high blood pressure) is defined as repeated measurement of a systolic value greater than 140 mmHg and/or a diastolic value greater than 90 mmHg.


Please note that this classification of blood pressure values is independent of age. Optimal blood pressure values have health benefits for all people.

There is no generally recognised definition of too-low blood pressure (hypotension). Readings of less than 100 mmHg systolic and less than 70 mmHg diastolic are considered too low. Please note that, unlike too-high blood pressure values, too-low blood pressure

values are not usually expected to be associated with health risks. However, if you are always feeling unwell, you should check with your doctor.

2.2 The importance of self-measurement of blood pressure

Constantly elevated blood pressure multiplies the risk for other health problems. Chief among these are the physical consequences of a heart attack or stroke including disability requiring nursing care, paralysis of one side of the body and organic heart or brain damage. In addition to other medical treatment measures, daily blood pressure monitoring will help to protect you from these adverse health outcomes.


2.3 Objectives of self-measurement of blood pressure

Your personal blood pressure profile is an important source of information. In case of drug treatment (e.g. in high blood pressure) your doctor can use your blood pressure profile to tailor your treatment to your particular needs. The better your drug treatment is tailored to your needs, the better you will feel, and the less you will suffer from side effects. Regular, accurate blood pressure

monitoring with Tensoval mobil will help you achieve this goal. Many people manage to lower their blood pressure through life-style changes (such as losing weight, dietary modification and getting more exercise) to levels that do not require drug treatment. Your blood pressure profile thus gives you excellent feedback on how successful your life-style changes are.

2.4 Regular blood pressure measurement


Numerous factors including physical exertion, taking drugs or the time of day may have an impact on blood pressure. Blood pressure should therefore always be measured at the same time of day under similar conditions. Our heart may beat up to 100,000 times a day, producing 100,000 different blood pressure values.


3. Getting ready for self-measurement

3.1 Inserting the batteries


Open the battery cover in the direction of the arrow on the upper side of the device by applying gentle pressure on the housing notches.

The marking in the battery compartment shows you which way to insert the batteries:

Looking for the “+” and “-” signs on the left and right ends of the batteries, insert the two “AAA” size batteries so that the positive (+) and negative (-) polarities match the “+” and “-” signs shown on the sticker in the battery compartment. If the polarities do not match, the device will not work and the batteries may leak!


Replace the lid of the battery compartment, applying gentle pressure until the hook snaps in place. The measured values recorded in the memory will not be lost when changing the batteries. The date and time must be reset.

3.2 Setting date and time

The batteries must have been inserted and the device switched off before you can set the date and time. Remove the cover of the batteries on the upper side of the device. Lift one of the inserted batteries for 10 seconds, and then push the batteries back into the compartment. Then you can close the battery compartment. Now the time function automatically appears. The display shows the number "31" for the day and "12" for the month. The date set is thus

the 31st of December.

The number on the left in the display flashes. Press the M1 (+) or M2 (-) buttons to change the day displayed. For example, press M2 (-) twice to set the date to the 29th of December. Store the current day by pressing the blue START/STOP button. Now the number on the right flashes. The display shows now the month. Also press the M1 (+) or M2 (-) buttons to display the current month and store it by pressing the blue START/STOP button. Now the year 2009 is displayed. You can change the year displayed as mentioned above and store it by pressing the START/STOP button (see fig.):


Then you may set the time. The number on the left in the display flashes, indicating 12:00. For ex-

ample, press M1 (+) twice to set the time to 14:00. Once the desired number of hours has been set, store it by pressing the START/STOP button. Now the number on the right flashes. Here you can change the minutes display. Once the desired number of minutes has been reached, store it by pressing the START/STOP button. All of your personal settings have now been saved completely.

3.3 Correct posture and position


- Measurements can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
 - The device must be placed at the level of the heart to give an accurate measurement result. To do this, hold the cuff at the level of the heart, supporting the arm at the elbow.
- Be sure to rest for 5 minutes before taking blood pressure measurements.
 - Do not talk or move during the measuring procedure. This may lead to false measurements.
 - You should take your blood pressure in a quiet place, in a relaxed seated position.
 - You should wait an hour before taking measurements after having coffee or smoking.
 - Do not take your blood pressure after taking a bath or exercising.
 - If you feel an urge to urinate, you should go to the toilette before taking measurements.

3.4 Applying the wrist monitor

Measurement should be taken on the naked wrist which usually has the higher blood pressure value. If you are unsure which arm to use,

ask your doctor. The device is permanently connected to the cuff, and the cuff should not be removed from the device. Now wrap the cuff round the wrist. The wrist monitor should be positioned on the inside of the wrist about 2 cm below the base of the wrist. The lettering on the upper side of the device points in your direction (see fig.):

The cuff should be firm but not too tight.

 Please note that inappropriate application of the cuff can lead to inaccurate readings. Also use the markings on the edge of the cuff to check if the cuff size is correct: The red dot on the silver-coloured strap should point on the red marking strip. If the red dot is outside the red marking strip, the


cuff is too small for you. The Tensoval mobil cuff circumference ranges from 12.5 to 22.5 cm. People with larger wrist circumferences are recommended to use an upper arm blood pressure monitor instead because wrist monitors give inaccurate blood pressure readings in such individuals (please read Chapter 7 “Tensoval range of products” to learn more about this).

4. Measuring blood pressure

We recommend carrying out blood pressure measurements whilst seated. Do not turn the device on until the cuff has been applied, otherwise the cuff can become damaged through the resulting excess pressure.

Press the blue START/STOP button. The appearance of all display segments followed by a flashing arrow pointing downwards, shows that the device is checking itself automatically and is ready for use.

Position the wrist with the strapped-on device and the palm of the hand pointing inwards at the level of the heart.

 **Important:** You should not move or talk throughout the entire measuring procedure! If you wish to stop measurement for any reason, simply press the blue START/STOP button. The inflation or measuring procedure is interrupted and an automatic fall in pressure occurs.

The Comfort Air Technology automatically determines the systolic blood pressure value and thus the individual inflation pressure for blood pressure measurement. If this inflation pressure is insufficient or if the measurement is interrupted, the device continues to pump at a rate of 40 mmHg until a high enough pressure is reached. Once the required inflation pressure has been reached, a short beep signals that the measurement starts.

If you require a higher inflation pressure, you can avoid having to repeat pumping by pressing the blue START/STOP button again shortly after inflation starts and holding down on it until the desired cuff pressure is reached. This should be approx. 30 mmHg over the systolic (upper) value.

As pressure in the cuff decreases, the heart symbol and the falling cuff pressure are displayed. A long beep indicates the end of measurement. Then the systolic and diastolic blood pressure values appear simultaneously on the display, with the pulse rate beneath them (see fig.):

Once the measurement is finished, the time appears above the measured values and M1 or M2 is displayed on the left. M1 represents

the measured values for a first person. Using M2 you can store the measured values for a second person (see 5.1 “Saving the measured values”).

In order to switch off the device, press the blue START/STOP button. Otherwise the device will switch itself off automatically after 3 minutes.

5. Setting the memory function

5.1 Saving the measured values

The device features two memory buttons, M1 and M2, to store readings for two different persons. M1 represents the measured values for a first person and M2 represents the measured values for a second person.

When a beep indicates the end of measurement, you may press M1 or M2 to apply the measured value to the respective person. This assignment can be made as long as the values are displayed. If you do not apply them, the measured value is automatically stored to the displayed memory.

5.2 Recalling the measured values

The device must be switched off to recall data from memory. Press M1 for the first person's stored values and M2 for the second person's values. The corresponding symbols M1 or M2 will be displayed. The corresponding symbols M1 or M2 will be displayed.

First, the mean value of all stored values for a person will be displayed. An “A” (for the designation “Average”) appears on the display and the number displayed in the top right-hand corner indicates the number of measurements, from which the average value was calculated (see fig.).

When pressing again, the most recently stored measured value will be displayed. By repeatedly pressing the memory button all stored values of the selected memory can be recalled one by one. The measured value and the number of the memory position are displayed when you access stored values. The number of the memory position, the date and time appear in 2 – 3 second intervals.

Tensoval mobil can store up to 60 measured values for each person (M1 or M2). The most recent measured value is always in memory position No.1. When all memory positions are occupied the oldest value will be deleted each time.

⚠ Important: The calculated mean value is based on all measured values stored for the respective person. If only two mea-

sured values are in memory, the mean value will be calculated from these two measured values. If only one measured value is in memory, no mean value will be calculated. You can cancel the memory's recall of data at any time by pressing the START/STOP button. Otherwise the device will switch itself off automatically after approx. 30 seconds. Even if the power supply fails, e.g. when changing the batteries, the stored values will still be available.

5.3 Deleting the measured values

You can delete all stored data either in M1 or in M2. To do this, press the memory button of the corresponding memories M1 or M2. The average value will appear on the display. Hold down the memory button for a longer time. After four seconds the display will flash and after further four seconds all data will be deleted. On the display appears only M1 or M2.

⚠ Important: If you release the memory button before the end of the time limit, no data will be deleted. If there is only a single


measured value in memory, this stored value cannot be deleted. It is not possible to delete individual values.

5.4 Using the guest mode

When a third person uses Tensoval mobil, it is recommended that they use the guest mode. This one is used to store a measured value neither to the memory positions M1 nor to M2. This avoids bias to the average values and stored values of the two main users of the device.

To take measurements in guest mode, start the measuring procedure by simultaneously pressing memory buttons M1 and M2. The START/STOP button needs not be pressed.

At the end of measurement the display simultaneously shows both M1 and M2 above the measured values rather than either symbol alone. The measured value can thus not be applied to either person, and the measured values will not be stored.

To switch off the device, press the START/STOP button, also in guest mode. Otherwise the device will switch itself off automatically after 3 minutes.

6. Maintenance of the device

Do not expose the device neither to extreme temperatures nor to humidity, dust or direct sunlight because this may lead to malfunction. This device consists of high-quality electronic precision components. Protect the device from shock and do not immerse in water. Only use a soft, moistened cloth to clean the device. Please use neither a diluter nor alcohol, detergents or solvents. The cuff can be cleaned carefully with a lightly moistened cloth and mild, pH neutral soap. Do not completely immerse the cuff in water and do not remove it from the device.

To protect them from external influences keep the wrist monitor as well as these instructions in the storage box.

7. Tensoval range of products

Apart from this wrist monitor, the HARTMANN range of products also includes upper arm blood pressure monitors. If you would like to learn more about other HARTMANN blood pressure devices, contact your specialist medical supplier (pharmacies or medical supply stores). Here is an overview of our product range:

- Upper arm blood pressure monitor Tensoval duo control featuring Duo Sensor Technology
- Upper arm blood pressure monitor Tensoval comfort
- Accessory cuffs such as:
 - Large cuff for upper arm circumference 32 – 42 cm
 - Standard draw-clamp cuff for upper arm circumference 22 – 32 cm
 - Preformed moulded cuff for upper arm circumference 22 – 32 cm

- HARTMANN mains adapter (suitable for upper arm blood pressure monitors only)

8. Warranty conditions

The device comes with a three-year warranty from the date of purchase. Accessories which are subject to wear and tear are not covered by this warranty (e.g. batteries, cuffs etc.). Further information about warranty conditions can be found in the separate warranty certificate included with this product.

9. Explanation of error displays

Error which has occurred	Possible causes	Remedy
Device will not turn on.	No batteries, they have been inserted incorrectly or are dead.	Check batteries and insert two identical, new batteries if necessary.
Cuff will not inflate.	Cuff defective.	Return the device to the designated service center for inspection.
	Moving or talking during the measuring procedure.	Do not talk or move during the measuring procedure.
	Cuff is too loosely applied.	Apply the cuff so that it fits snugly around the wrist.
	Cuff loses air/has a leak.	Return the device to the designated service address.
	Measurement inaccurate.	Please rest for at least a minute and take the measurement again.
	The pressure in cuff exceeds 300 mmHg. Therefore, an automatic fall in pressure occurs.	Please rest for at least a minute and take the measurement again.

Error which has occurred	Possible causes	Remedy
	Batteries are almost dead. Only a few more measurements are possible (approx. 30).	Keep new, identical batteries handy (type "AAA" LR03).
	Batteries are dead and have to be replaced.	Insert new, identical batteries (type "AAA" LR03).
Implausible measured values.	Device not positioned at the level of the heart.	Position the wrist at the level of the heart, and repeat the measurement.
	Incorrect cuff size.	Use an upper arm blood pressure monitor.
	Cuff placed on top of clothing.	Apply cuff on the naked skin.
	Rolled-up clothing impedes blood circulation.	Wear loose clothing. Make sure that rolled-up sleeves do not impair circulation in the arm.
	Wrist jewellery and/or a watch prevent(s) blood from circulating.	Take off your watch and/or any wrist jewellery you may be wearing before measuring your blood pressure.


Error which has occurred	Possible causes	Remedy
Implausible measured values.	Moving, talking or excitement during the measuring procedure.	Please take measurements in a relaxed position whilst seated. Do not talk or move during the measuring procedure.
	Lack of relaxation during taking a measurement.	Relax for about 5 minutes before taking a measurement.
	Stimulants taken before measurement.	Please avoid alcohol/ nicotine and caffeine for one hour before taking a measurement.
Memory cannot be deleted.	Only one stored value is in memory. The memory cannot be deleted until two values have been stored.	Save a second value. Then delete the memory.
	The memory can only be deleted when the average value is displayed.	Press the memory button until the average value is displayed.

Please check the possible causes, when an error symbol appears and note the instructions on self-measurement in Chapter 4. Relax for a minute and take the measurement again.

10. Important notes

10.1 Drugs

Self-measurement of blood pressure does not replace treatment! So do not interpret your measured values on your own and do not use them for self-prescribed treatment. Take measurements as instructed by your doctor and have confidence in your doctor's diagnosis. Take drugs as prescribed by your doctor and never alter the dose on your own. Discuss the appropriate time for self-measurement of blood pressure with your doctor.

10.2 Pregnancy

Blood pressure may change during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because the elevated blood pressure may affect the development of the foetus. Check with your doctor whether and, if so, when you should carry out self-measurement of blood pressure.

10.3 Diabetes, history of other medical conditions

If you have diabetes or narrowed blood vessels (arteriosclerosis) you should consult your doctor before carrying out self-measurement because greatly deviating individual measured values may occur in such cases.

10.4. Arrhythmias, heart rhythm disorders (cardiac arrhythmias), cardiac pacemakers

Heart rhythm disorders (arrhythmias) are disturbances of the normal rate or rhythm of the heart-beat. A distinction should be made between mild and severe heart rhythm disorders. This can only be determined by a special investiga-

tion carried out by a doctor. If you have a heart rhythm disorder, you are specifically recommended to use our innovative upper arm blood pressure monitor Tensoval duo control. Using Duo Sensor Technology – a technology based on Korotkoff's principle – this device can sense various types of cardiac arrhythmias and give accurate readings. Please bear in mind: Severe heart rhythm disorders may produce false measurements or impair the measuring precision. Please discuss with your doctor whether self-measurement of blood pressure – and if yes – which type is suitable for you. Self-measurement of blood pressure may fail in cardiac pacemaker wearers, but the blood pressure device itself has no impact on the cardiac pacemaker. Please note that the displayed pulse rate is not suitable for checking the rate of cardiac pacemakers. Please check with your doctor whether self-measurement of blood pressure is advisable if you are wearing a cardiac pacemaker.

11. Control displays and symbols

Flashes when the device is measuring and the pulse is being taken.

Electronic device disposal note

PAUL HARTMANN AG
Marketer of this medical device

Change batteries!

Measuring error, cf. Chap. 9

Disposal note

Appears during inflation.

Appears during automatic checking.

Displays the stored measured values for person 1.

Displays the stored measured values for person 2.

Protection from electric shock (type BF).

Pay attention to the operating instructions.

12. Technical data

Measuring method:	Oscillometric
Display range:	0 – 297 mmHg
Measuring range:	Systole (SYS): 50 – 250 mmHg Diastole (DIA): 40 – 180 mmHg Pulse: 40–160 beats / minute
Technical measuring precision:	Cuff pressure: ± 3 mmHg, Pulse: ± 5 % of displayed pulse rate
Power supply:	2 x 1.5 V Mignon alkaline-manganese (AAA/LR03) batteries
Battery capacity:	1,000 measurements
Inflation pressure:	At least 150 mmHg
Automatic switch-off function:	3 minutes after end of measurement
Cuff:	12.5 – 22.5 cm
Pressure release valve:	Electronically controlled linear valve
Memory capacity:	2 x 60 measurements and mean value
Operating conditions:	Ambient temperature: $+10$ °C to $+40$ °C ($+ 50$ °F to $+ 104$ °F) Relative humidity: 15 – 90 %
Storage / transport conditions:	Ambient temperature: -20 °C to $+50$ °C ($- 4$ °F to $+122$ °F) Relative humidity: 15 – 90 %
Serial number:	See battery compartment

13. Power supply, disposal note and safety information

13.1 Batteries and disposal

- The two high-quality batteries included with the device guarantee about 1,000 measurements. Use only high-quality batteries (see specification in Chapter 12 “Technical data”). If you are using less efficient batteries we can no longer guarantee 1,000 measurements.
- Never mix old and new batteries or batteries made by different manufacturers.
- Immediately remove exhausted batteries.
- If you do not intend to use the device for a longer period of time, you should remove the batteries to prevent possible leakage.
- Please help protect the environment: Do not dispose of batteries in household waste. Use designated collection points or municipal collection/recycling centers to dispose of hazardous household waste.

- Information concerning the disposal of electronic devices (private households):

This symbol on products and/or accompanying documents means that spent electronic products must not be mixed with regular household waste. Take these products to designated collection points where they are accepted free of charge for proper treatment, reclamation and recycling. Proper disposal is in the interests of environmental protection and prevents possible harmful effects on people and the environment that may result from inappropriate handling of wastes. More detailed information about your nearest collection point is available from your local council office.

13.2 Safety information

- Do not leave the device unattended near toddlers or persons who cannot operate it themselves.
- Use the device for taking blood pressure measurements on the wrist only.

- Do not expose the device to hard knocks or vibrations.
- Do not drop the device to the floor. Do not excessively bend or fold the cuff.
- The device must not be altered, dismantled or repaired by the user.

14. Legal requirements and guidelines

Tensoval mobil complies with the requirements of the EC directive 93/42/EEC on medical devices (Medical Device Directive MDD) and bears the CE mark.

The device complies, for example, with the European Standard EN 1060 relating to non-invasive blood pressure measuring devices, Part 1: General requirements and Part 3: Additional requirements for electro-mechanical blood pressure measuring systems. Clinical testing of measurement precision was performed according to the European Standard EN 1060-4.

Manufacturer: PAUL HARTMANN AG,
D-89522 Heidenheim, Germany.

15. Measurement function check and service addresses

15.1 Explanation of measurement function check

Each Tensoval mobil device has been carefully tested by HARTMANN for measuring precision, and been developed for a long useable service life.

We recommend a check of the measurement function at intervals of two years in the case of professional users, such as pharmacies, medical practices or hospitals. You should also observe the national regulations, such as, in Germany, the "Medizinprodukte-Betreiberverordnung" (Medical Device Operating Regulation). Measurement function checks can be carried out either by the indicated Service address, by competent authorities or authorised maintenance providers against compensation

Instructions for the measurement function check:

A device function check can be carried out on people or using a suitable simulator. Measurement function check involves testing for

leak tightness of the pressure system and possible deviations of the pressure reading.
Remove at least one battery in order to switch to calibration mode. Hold down on the START/STOP button and then insert the battery. Release the button and, after a few moments, two zeros will appear one above each other on the display.
Instructions on the measurement function check will be supplied on request to competent authorities or authorised maintenance providers by HARTMANN.

15.2 Contact information for customer queries

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